


Anti-Alox5 Rabbit Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P09917

Item Number: CM0023979



Introduction

High-quality anti-ALOX5 rabbit polyclonal antibody validated for WB, IF/ICC, IHC-P, and ELISA. Recognizes human, mouse, and rat 5-lipoxygenase (P09917), a key enzyme in leukotriene biosynthesis and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	ALOX5 Rabbit pAb
Remarks/Alias	5-LO; 5LPG; LOG5; 5-LOX; ALOX5
Species	Human
Gene ID (Human)	240
Gene ID	240
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of ALOX5 (NP_000689.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09917
Protein Weight	78 kDa
Shipping	Ice bag