

Anti-Aloxe3 Polyclonal Antibody For Wb, Elisa - Q9Byj1

Item Number: CM0027075



Introduction

High-quality rabbit polyclonal antibody against human ALOXE3 (Epidermis-type lipoxygenase 3), validated for WB and ELISA with cross-reactivity to mouse and rat. Ideal for skin barrier and lipid metabolism research.

[Learn More](#)

Parameter	Value
Product Name	ALOXE3 Rabbit pAb
Remarks/Alias	LI5; ARC13; E-LOX; eLOX3; eLOX-3; ALOXE3
Species	Human
GeneID (Human)	59344
GeneID	59344
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human ALOXE3 (NP_067641.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BYJ1
Protein Weight	81kDa
Shipping	Ice bag