

Anti-Alpha-1 Antitrypsin (Serpina1) Monoclonal Antibody For Wb, If/Icc, Elisa - P01009

Item Number: CM0005644



Introduction

Rabbit monoclonal antibody against human Alpha-1 Antitrypsin (SERPINA1). Recombinant immunogen (aa 25-418). Validated for WB, IF/ICC, ELISA. Use in respiratory disease and serpin research.

[Learn More](#)

Parameter	Value
Product Name	Alpha-1 Antitrypsin (SERPINA1) Rabbit mAb
Remarks/Alias	PI; A1A; AAT; PI1; A1AT; nNIF; PRO2275; alpha1AT; Alpha-1 Antitrypsin (SERPINA1)
Species	Human
Gene ID (Human)	5265
Gene ID	5265
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-418 of human Alpha-1 Antitrypsin (SERPINA1) (NP_000286.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P01009
Protein Weight	47 kDa
Shipping	Ice bag