

Anti-Alpha-1 Antitrypsin (Serpina1) Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P01009

Item Number: CM0017339



Introduction

Alpha-1 Antitrypsin (SERPINA1) Rabbit pAb is a polyclonal antibody detecting human and mouse SERPINA1. Validated for WB, IHC-P, IF/ICC, and ELISA, it targets a key serine protease inhibitor involved in lung protection. Ideal for respiratory and inflammation research.

[Learn More](#)

Parameter	Specification
Product Name	Alpha-1 Antitrypsin (SERPINA1) Rabbit pAb
Remarks/Alias	PI; A1A; AAT; PI1; A1AT; nNIF; PRO2275; alpha1AT; Alpha-1 Antitrypsin (SERPINA1)
Species	Human
GeneID (Human)	5265
GeneID	5265
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 156-255 of human Alpha-1 Antitrypsin (SERPINA1) (NP_000286.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P01009
Protein Weight	47kDa
Shipping	Ice bag