


Anti-Alpha 1 Spectrin Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - P02549

Item Number: CM0017585



Introduction

Rabbit polyclonal antibody targeting human alpha 1 Spectrin (SPTA1), validated for WB, IF-P, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for cytoskeleton and erythrocyte membrane research.

[Learn More](#)

Parameter	Value
Product Name	alpha 1 Spectrin Rabbit pAb
Remarks/Alias	EL2; HPP; HS3; SPH3; SPTA; alpha 1 Spectrin
Species	Human
GeneID (Human)	6708
GeneID	6708
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 940-1160 of human alpha 1 Spectrin (NP_003117.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P02549
Protein Weight	280kDa
Shipping	Ice bag