


Anti-Alpha 1 Spectrin Rabbit Monoclonal Antibody For Wb & Elisa - P02549

Item Number: CM0010988



Introduction

Rabbit monoclonal antibody against human alpha 1 Spectrin (SPTA1, P02549). Validated for Western blot and ELISA, with cross-reactivity to mouse and rat. Essential for erythrocyte membrane cytoskeleton research.

[Learn More](#)

Parameters	Specifications
Product Name	alpha 1 Spectrin Rabbit mAb
Remarks/Alias	EL2; HPP; HS3; SPH3; SPTA; alpha 1 Spectrin
Species	Human
GeneID (Human)	6708
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 2300-2400 of human alpha 1 Spectrin (P02549).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P02549
Protein Weight	280kDa
Shipping	Ice bag