

Anti-Alpha-2-Macroglobulin (A2M) Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P01023

Item Number: CM0008560



Introduction

High-quality rabbit monoclonal antibody targeting human Alpha-2-Macroglobulin (A2M), a key protease inhibitor. Validated for western blot, immunofluorescence, immunohistochemistry paraffin, and ELISA. Works across human, mouse, and rat species. Ideal for studying proteinase inhibition mechanisms.

[Learn More](#)

Parameter	Value
Product Name	Alpha-2-Macroglobulin (A2M) Rabbit mAb
Remarks/Alias	A2MD, CPAMD5, FWP007, S863-7, Alpha-2-Macroglobulin (A2M)
Species	Human
GeneID (Human)	2
GeneID	2
Immunogen	Synthetic peptide corresponding to amino acids 1200-1300 of human A2M (P01023)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01023
Protein Weight	163 kDa
Shipping	Ice bag