

Anti-Alpha-2-Macroglobulin (A2M) Polyclonal Antibody For Wb - P01023

Item Number: CM0023080



Introduction

Rabbit polyclonal antibody against human Alpha-2-Macroglobulin (A2M), validated for WB, IHC-P, and ELISA. Useful for investigating protease inhibition mechanisms and plasma protein interactions.

[Learn More](#)

Parameter	Value
Product Name	Alpha-2-Macroglobulin (A2M) Rabbit pAb
Remarks/Alias	A2MD; CPAMD5; FWP007; S863-7; Alpha-2-Macroglobulin (A2M)
Species	Human
GeneID	2
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1155-1474 of human Alpha-2-Macroglobulin (A2M) (NP_000005.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P01023
Protein Weight	163kDa
Shipping	Ice bag