


## Anti-Alpha + Beta Synuclein Monoclonal Antibody For Wb, If-P, Elisa - P37840 / Q16143

Item Number: CM0010481



### Introduction

Rabbit monoclonal antibody against human alpha + beta Synuclein (P37840/Q16143) for Western blot, immunofluorescence, ELISA. Cross-reacts with mouse, rat. Ideal for Parkinson's and neurodegeneration research.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | alpha + beta Synuclein Rabbit mAb                                                                                  |
| Remarks/Alias    | NACP; PARK1; PARK4; PD1; alpha + beta Synuclein                                                                    |
| Species          | Human                                                                                                              |
| GeneID (Human)   | 6620/6622                                                                                                          |
| GeneID           | 6622/6620                                                                                                          |
| Immunogen        | A synthetic peptide corresponding to a sequence within amino acids 1-100 of human alpha + beta Synuclein (P37840). |
| Source           | Rabbit                                                                                                             |
| Category         | Monoclonal Antibodies                                                                                              |
| Application      | WB, IF-P, ELISA                                                                                                    |
| Cross Reactivity | Human, Mouse, Rat                                                                                                  |
| SWISS            | P37840/Q16143                                                                                                      |
| Protein Weight   | 18kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                            |