

Anti-Alpha-Fetoprotein (Afp) Monoclonal Antibody For Wb, Elisa - P02771

Item Number: CM0006798



Introduction

Recombinant rabbit monoclonal antibody specific for human Alpha-Fetoprotein (AFP) protein. Validated for Western blot (WB) and ELISA applications. Detects a ~69 kDa band. Ideal for cancer biomarker and prenatal diagnostic research.

[Learn More](#)

Parameter	Value
Product Name	Alpha-Fetoprotein (AFP) Rabbit mAb
Remarks/Alias	AFPD; FETA; HPAFP; Alpha-Fetoprotein (AFP)
Species	Human
GeneID (Human)	174
GeneID	174
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P02771
Protein Weight	69kDa
Shipping	Ice bag