

Anti-Alpha-Fetoprotein (Afp) Monoclonal Antibody For Wb, If/Icc, Elisa - P02771

Item Number: CM0004733



Introduction

Recombinant rabbit monoclonal antibody against human Alpha-Fetoprotein (AFP), validated for WB, IF/ICC, and ELISA. Ideal for liver cancer biomarker research and diagnostic assay development.

[Learn More](#)

Parameter	Specification
Product Name	Alpha-Fetoprotein (AFP) Rabbit mAb
Remarks/Alias	AFPD; FETA; HPAFP; Alpha-Fetoprotein (AFP)
Species	Human
GeneID (Human)	174
GeneID	174
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-609 of human Alpha-Fetoprotein (AFP) (NP_001125.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P02771
Protein Weight	69kDa
Shipping	Ice bag