

Anti-Alpha-Fetoprotein (Afp) Monoclonal Antibody For Wb, If/Icc, Elisa - P02771

Item Number: CM0003326



Introduction

Mouse monoclonal antibody targeting human Alpha-Fetoprotein (AFP). Suitable for Western blot, IF/ICC, and ELISA applications. Cross-reacts with human, mouse, and rat samples. Ideal for liver cancer biomarker research.

[Learn More](#)

Parameter	Value
Product Name	Alpha-Fetoprotein (AFP) Mouse mAb
Remarks/Alias	AFPD; FETA; HPAFP; Alpha-Fetoprotein (AFP)
Species	Human
GeneID (Human)	174
GeneID	174
Immunogen	Recombinant protein of human Alpha-Fetoprotein (Alpha-Fetoprotein (AFP))
Source	Mouse
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P02771
Protein Weight	69kDa
Shipping	Ice bag