

Anti-Alpha-Fetoprotein (Afp) Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P02771

Item Number: CM0014059



Introduction

Rabbit polyclonal antibody against human Alpha-Fetoprotein (AFP), suitable for WB, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Detects a 69 kDa band. Ideal for developmental biology and oncology research.

[Learn More](#)

Parameter	Value
Product Name	Alpha-Fetoprotein (AFP) Rabbit pAb
Remarks/Alias	AFPD; FETA; HPAFP; Alpha-Fetoprotein (AFP)
Species	Human
GeneID (Human)	174
GeneID	174
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 360-609 of human Alpha-Fetoprotein (AFP) (NP_001125.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P02771
Protein Weight	69kDa
Shipping	Ice bag