


Anti-Alpha Fodrin Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q13813

Item Number: CM0017102



Introduction

Rabbit polyclonal antibody targeting human Alpha Fodrin (SPTAN1). Confirmed for WB, IHC-P, IF/ICC, and ELISA; cross-reacts with mouse and rat. Ideal for cytoskeleton and secretion research.

[Learn More](#)

Parameter	Value
Product Name	Alpha Fodrin Rabbit pAb
Remarks/Alias	DEE5; NEAS; EIEE5; SPTA2; Alpha Fodrin
Species	Human
GeneID (Human)	6709
GeneID	6709
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 950-1130 of human Alpha Fodrin (NP_001123910.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13813
Protein Weight	285kDa
Shipping	Ice bag