


Anti-Alpha Fodrin Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q13813

Item Number: CM0007604



Introduction

High-quality rabbit monoclonal antibody against human Alpha Fodrin (SPTAN1). Validated for WB, IHC-P, ELISA with cross-reactivity to mouse and rat. Ideal for cytoskeletal research and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	Alpha Fodrin Rabbit mAb
Remarks/Alias	DEE5; NEAS; EIEE5; SPTA2; Alpha Fodrin
Species	Human
Gene ID (Human)	6709
Gene ID	6709
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Alpha Fodrin (Q13813).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13813
Protein Weight	285kDa
Shipping	Ice bag