

Anti-Alpha Internexin Polyclonal Antibody For Wb, Elisa - Q16352

Item Number: CM0027595

Introduction

Alpha Internexin (NF-66, NEF5) rabbit polyclonal antibody validated for WB and ELISA. Recognizes human and mouse alpha-internexin, a neuronal intermediate filament protein. Suitable for neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	Alpha Internexin Rabbit pAb
Remarks/Alias	NEF5; NF66; NF-66; TXBP-1; Alpha Internexin
Species	Human
GeneID (Human)	9118
GeneID	9118
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 13-300 of human Alpha Internexin (Q16352)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q16352
Protein Weight	55kDa
Shipping	Ice bag