

Anti-Alpha Internexin Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q16352

Item Number: CM0008489

Introduction

High-quality rabbit monoclonal antibody against human Alpha Internexin (INA, Q16352), validated for WB, IF-P, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for neuronal intermediate filament research.

[Learn More](#)

Parameter	Value
Product Name	Alpha Internexin Rabbit mAb
Remarks/Alias	NEF5; NF66; NF-66; TXBP-1; Alpha Internexin
Species	Human
Gene ID (Human)	9118
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 350-450 of human Alpha Internexin (Q16352).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16352
Protein Weight	55kDa
Shipping	Ice bag