

Anti-Alpha Sarcoglycan (Sgca) Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q16586

Item Number: CM0003401

Introduction

Rabbit monoclonal antibody against human alpha-sarcoglycan (SGCA, 50DAG, adhalin). Validated for WB, IF-P, IHC-P, ELISA; cross-reacts with mouse, rat. Targets a 43 kDa dystrophin-associated glycoprotein critical for muscle membrane stability.

[Learn More](#)

Parameter	Value
Product Name	alpha Sarcoglycan (SGCA) Rabbit mAb
Remarks/Alias	ADL; DAG2; 50DAG; DMDA2; LGMD2D; LGMDR3; SCARMD1; adhalin; alpha Sarcoglycan (SGCA)
Species	Human
GeneID (Human)	6442
GeneID	6442
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q16586
Protein Weight	43kDa
Shipping	Ice bag