

Anti-Alpha Sarcoglycan (Sgca) Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q16586

Item Number: CM0026439

Introduction

High-quality rabbit polyclonal antibody against human alpha-sarcoglycan (SGCA), validated for WB, IF/ICC, and ELISA. Cross-reacts with human, mouse, rat. Ideal for muscle disease and sarcolemma research.

[Learn More](#)

Parameter	Value
Product Name	alpha Sarcoglycan (SGCA) Rabbit pAb
Remarks/Alias	ADL; DAG2; 50DAG; DMDA2; LGMD2D; LGMDR3; SCARMD1; adhalin; alpha Sarcoglycan (SGCA)
Species	Human
GeneID (Human)	6442
GeneID	6442
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 60-290 of human alpha Sarcoglycan (alpha Sarcoglycan (SGCA)) (NP_000014.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16586
Protein Weight	43kDa
Shipping	Ice bag