

## Anti-Alpha Snap Monoclonal Antibody For Wb And Elisa - P54920

Item Number: CM0003174

### Introduction

Recombinant rabbit monoclonal antibody targeting human Alpha SNAP (NAPA). Validated for WB and ELISA; cross-reacts with mouse and rat. Suitable for vesicular transport and membrane trafficking research.

[Learn More](#)

| Parameter        | Specification                                                                                                      |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | Alpha SNAP Rabbit mAb                                                                                              |
| Remarks/Alias    | SNAPA; Alpha SNAP                                                                                                  |
| Species          | Human                                                                                                              |
| Gene ID (Human)  | 8775                                                                                                               |
| Gene ID          | 8775                                                                                                               |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 201-295 of human Alpha SNAP (P54920) |
| Source           | Rabbit                                                                                                             |
| Category         | Monoclonal Antibodies                                                                                              |
| Application      | WB, ELISA                                                                                                          |
| Cross Reactivity | Human, Mouse, Rat                                                                                                  |
| SWISS            | P54920                                                                                                             |
| Protein Weight   | 33kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                            |