

Anti-Alpha Snap Polyclonal Antibody For Wb, Ihc-P, Elisa - P54920

Item Number: CM0030604

Introduction

CamelBio's Alpha SNAP Rabbit Polyclonal Antibody (P54920) is intended for WB, IHC-P, and ELISA. It cross-reacts with human, mouse, and rat Alpha SNAP, a key mediator of ER-Golgi vesicular transport.

[Learn More](#)

Parameter	Value
Product Name	Alpha SNAP Rabbit pAb
Remarks/Alias	SNAPA; Alpha SNAP
Species	Human
GeneID (Human)	8775
GeneID	8775
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-295 of human Alpha SNAP (NP_003818.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P54920
Protein Weight	33kDa
Shipping	Ice bag