

Anti-Alpha-Synuclein Monoclonal Antibody For Wb, If, Ihc-P, Elisa - P37840

Item Number: CM0002638



Introduction

This rabbit monoclonal antibody against human alpha-synuclein (SNCA). Validated for WB, IF, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for neuroscience and neurodegeneration research.

[Learn More](#)

Parameter	Value
Product Name	α -Synuclein Rabbit PolymAb®
Remarks/Alias	PD1; NACP; PARK1; PARK4
Species	Human
Gene ID	6622
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 59-140 of human α -Synuclein (NP_000336.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-F, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P37840
Protein Weight	14 kDa
Shipping	Ice bag