


## Anti-Alpha-Synuclein Truncation 1-119 Monoclonal Antibody [Mjff-C-24D2] For Slot Blot, If/Icc, Elisa - P37840

Item Number: CM0008622



### Introduction

Rabbit monoclonal antibody (clone MJFF-C-24D2) targeting human and mouse alpha-synuclein truncation 1-119. Validated in Slot blot, IF/ICC, and ELISA. Suitable for studying Parkinson's disease and protein aggregation in neuroscience research.

[Learn More](#)

| Parameter        | Value                                                                                                                      |
|------------------|----------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Alpha-synuclein Truncation 1-119 mAb [MJFF-C-24D2]                                                                         |
| Remarks/Alias    | PD1; NACP; PARK1; PARK4                                                                                                    |
| Species          | Human, Mouse                                                                                                               |
| GeneID (Human)   | 6622                                                                                                                       |
| GeneID           | 6622                                                                                                                       |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-119 of human alpha-synuclein (NP_000336.1) |
| Source           | Rabbit                                                                                                                     |
| Category         | Monoclonal Antibodies                                                                                                      |
| Application      | Slot blot, IF/ICC, ELISA                                                                                                   |
| Cross Reactivity | Human, Mouse                                                                                                               |
| SWISS            | P37840                                                                                                                     |
| Protein Weight   | 14kDa                                                                                                                      |
| Shipping         | Ice bag                                                                                                                    |