


Anti-Alpha-Synuclein Truncation 1-122 Monoclonal Antibody For Slot Blot, Elisa - P37840

Item Number: CM0006210



Introduction

Rabbit monoclonal antibody targeting human alpha-synuclein truncation 1-122. Validated for slot blot and ELISA, cross-reacts with human. Useful for studying Parkinson's disease, synaptic function, and amyloid pathology.

[Learn More](#)

Parameter	Value
Product Name	Alpha-synuclein Truncation 1-122 mAb [MJFF-D-10B7]
Remarks/Alias	PD1; NACP; PARK1; PARK4
Species	Human
GeneID (Human)	6622
GeneID	6622
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-122 of human alpha-synuclein (NP_000336.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	Slot blot, ELISA
Cross Reactivity	Human
SWISS	P37840
Protein Weight	14kDa
Shipping	Ice bag