

Anti-Alpi (Iap) Polyclonal Antibody For Wb, Elisa - P09923

Item Number: CM0028909

Introduction

ALPI Rabbit Polyclonal Antibody validated for WB and ELISA, cross-reactive with human and mouse. Targets intestinal alkaline phosphatase (IAP), a GPI-anchored membrane hydrolase involved in phosphate metabolism and gut homeostasis.

[Learn More](#)

Parameter	Value
Product Name	ALPI Rabbit pAb
Remarks/Alias	IAP; ALPI
Species	Human
Gene ID	248
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P09923
Protein Weight	57kDa
Shipping	Ice bag