

Anti-Alpk1 Rabbit Polyclonal Antibody For Wb And Elisa - Q96Qp1

Item Number: CM0019498



Introduction

High-quality rabbit polyclonal antibody targeting human ALPK1. Validated for Western blot and ELISA, this antibody detects the innate immune kinase ALPK1, which senses bacterial ADP-heptose and activates NF-κB signaling.

[Learn More](#)

Parameter	Value
Product Name	ALPK1 Rabbit pAb
Remarks/Alias	LAK; ROSAH; 8430410J10Rik; ALPK1
Species	Human
Gene ID	80216
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 340-500 of human ALPK1 (NP_001095876.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q96QP1
Protein Weight	139kDa
Shipping	Ice bag