

Anti-ALPK1 Rabbit Polyclonal Antibody For Wb - Q96Qp1

Item Number: CM0019643



Introduction

High-quality rabbit polyclonal antibody targeting human ALPK1, validated for Western blot. Recognizes alpha-protein kinase 1, a key innate immune kinase sensing bacterial ADP-Heptose and activating NF- κ B. Ideal for infection and inflammation studies.

[Learn More](#)

Parameter	Value
Product Name	ALPK1 Rabbit pAb
Remarks/Alias	LAK; ROSAH; 8430410J10Rik; ALPK1
Species	Human
Gene ID	80216
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-473 of human ALPK1 (NP_079420.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB
Cross Reactivity	Human
SWISS	Q96QP1
Protein Weight	139kDa
Shipping	Ice bag