

Anti-ALPPL2 Rabbit Polyclonal Antibody For Wb, Elisa - P10696

Item Number: CM0025743

Introduction

Rabbit polyclonal antibody against human ALPPL2 (GCAP) for WB and ELISA. Cross-reacts with human and mouse. Ideal for germ cell tumor research and IVD assay development. Protein weight ~57 kDa.

[Learn More](#)

Parameter	Value
Product Name	ALPPL2 Rabbit pAb
Remarks/Alias	GCAP, ALPPL, ALPPL2
Species	Human
GeneID (Human)	251
GeneID	251
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-280 of human ALPPL2 (NP_112603.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P10696
Protein Weight	57kDa
Shipping	Ice bag