

Anti-Als2Cr1 Rabbit Monoclonal Antibody For Wb, Elisa - Q9Gzt8

Item Number: CM0008900

Introduction

High-quality rabbit monoclonal antibody targeting human ALS2CR1 (NIF3-like protein 1), validated for Western blot and ELISA. Suitable for neuronal differentiation and gene expression studies.

[Learn More](#)

Parameter	Value
Product Name	ALS2CR1 Rabbit mAb
Remarks/Alias	CALS-7; MDS015; ALS2CR1
Species	Human
GeneID (Human)	60491
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9GZT8
Protein Weight	42 kDa
Shipping	Ice bag