


Anti-Alx4 Rabbit Polyclonal Antibody For Wb, Elisa - Q9H161

Item Number: CM0030920



Introduction

High-quality rabbit polyclonal antibody against ALX4, a homeobox transcription factor involved in craniofacial and limb development. Validated for WB and ELISA, with cross-reactivity to human, mouse, and rat. Ideal for developmental biology research.

[Learn More](#)

Parameter	Value
Product Name	ALX4 Rabbit pAb
Remarks/Alias	CRS5; FND2; ALX4
Species	Human
GeneID (human)	60529
GeneID	60529
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human ALX4 (NP_068745.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H161
Protein Weight	44kDa
Shipping	Ice bag