

Anti-Amacr Rabbit Polyclonal Antibody For Wb, If, Ihc, Elisa - Q9Uhk6

Item Number: CM0017990

Introduction

Rabbit polyclonal antibody targeting human AMACR (Alpha-methylacyl-CoA racemase), suited for Western blot, immunofluorescence, immunohistochemistry, and ELISA applications. Cross-reacts with mouse and rat. Valuable tool for investigating peroxisomal fatty acid metabolism, bile acid synthesis, and related metabolic disorders.

[Learn More](#)

Parameter	Value
Product Name	AMACR Rabbit pAb
Remarks/Alias	RM; RACE; CBAS4; P504S; AMACRD; AMACR
Species	Human
GeneID (Human)	23600
GeneID	23600
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-138 of human AMACR (NP_055139.4).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UHK6
Protein Weight	42kDa
Shipping	Ice bag