


Anti-Ambp Rabbit Polyclonal Antibody For Wb, Elisa - P02760

Item Number: CM0024164



Introduction

High-quality rabbit polyclonal antibody targeting human AMBP protein. Suitable for Western blot (WB) and ELISA applications; cross-reacts with mouse. Ideal for studying oxidative stress, protease inhibition, and tissue repair processes.

[Learn More](#)

Parameter	Value
Product Name	AMBP Rabbit pAb
Remarks/Alias	A1M, HCP, ITI, UTI, EDC1, HI30, ITIL, IATIL, ITILC, AMBP
Species	Human
GeneID (human)	259
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-202 of human AMBP (NP_001624.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P02760
Protein Weight	39kDa
Shipping	Ice bag