


Anti-Ambra1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9C0C7

Item Number: CM0017093



Introduction

Rabbit polyclonal antibody against human Ambra1, suitable for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Ideal for autophagy, cell cycle, and tumor suppression research.

[Learn More](#)

Parameter	Value
Product Name	Ambra1 Rabbit pAb
Remarks/Alias	DCAF3; WDR94; Ambra1
Species	Human
Gene ID (Human)	55626
Gene ID	55626
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1109-1208 of human Ambra1 (NP_060219.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9C0C7
Protein Weight	143 kDa
Shipping	Ice bag