

Anti-Amd1 Polyclonal Antibody For Wb, Elisa - P17707

Item Number: CM0028469



Introduction

AMD1 Rabbit pAb specifically detects S-adenosylmethionine decarboxylase proenzyme (SAMDC/ADOMETDC) in human, mouse, and rat samples by Western blot and ELISA. Ideal for studying polyamine metabolism and embryonic stem cell maintenance.

[Learn More](#)

Parameter	Value
Product Name	AMD1 Rabbit pAb
Remarks/Alias	AMD; SAMDC; ADOMETDC; AMD1
Species	Human
GeneID (Human)	262
GeneID	262
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 75-334 of human AMD1 (NP_001625.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17707
Protein Weight	38kDa
Shipping	Ice bag