


Anti-Amfr Rabbit Polyclonal Antibody For Wb And Elisa - Q9Ukv5

Item Number: CM0026568



Introduction

Rabbit polyclonal antibody targeting human AMFR (gp78), an E3 ubiquitin-protein ligase involved in ER-associated degradation. Validated for WB and ELISA. Ideal for ubiquitination and ER stress research.

[Learn More](#)

Parameter	Value
Product Name	AMFR Rabbit pAb
Remarks/Alias	GP78; RNF45; AMFR
Species	Human
Gene ID (Human)	267
Gene ID	267
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 384-643 of human AMFR (NP_001135.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9UKV5
Protein Weight	73kDa
Shipping	Ice bag