


Anti-Amh Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - P03971

Item Number: CM0024727



Introduction

Rabbit polyclonal antibody against human AMH (Muellerian-inhibiting factor/MIS). Validated for WB, IF-P, IHC-P, ELISA, and reacts with human, mouse, rat. Ideal for reproductive biology and developmental studies.

[Learn More](#)

Parameter	Value
Product Name	AMH Rabbit pAb
Remarks/Alias	MIF; MIS; AMH
Species	Human
GeneID (Human)	268
GeneID	268
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 360-560 of human AMH (NP_000470.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P03971
Protein Weight	59kDa
Shipping	Ice bag