


Anti-Amhr2 Monoclonal Antibody For Wb, Elisa - Q16671

Item Number: CM0007799



Introduction

Rabbit monoclonal antibody against human AMHR2 (Anti-Muellerian hormone type-2 receptor), aliases AMHR/MRII/MISR2/MISRII. Validated for WB and ELISA. Immunogen: synthetic peptide within aa 50-150. For reproductive biology and signaling research.

[Learn More](#)

Parameter	Value
Product Name	AMHR2 Rabbit mAb
Remarks/Alias	AMHR; MRII; MISR2; MISRII; AMHR2
Species	Human
GeneID (Human)	269
GeneID	269
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 50-150 of human AMHR2 (NP_065434.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q16671
Protein Weight	52kDa/63kDa
Shipping	Ice bag