


Anti-Amhr2 Rabbit Monoclonal Antibody For Wb, Elisa - Q16671

Item Number: CM0006529



Introduction

Rabbit monoclonal antibody against human AMHR2 (Anti-Muellerian hormone type-2 receptor). Validated for WB and ELISA. Ideal for studying TGF-beta signaling, sex determination, and reproductive biology. Recognizes human AMHR2 protein.

[Learn More](#)

Parameter	Value
Product Name	AMHR2 Rabbit mAb
Remarks/Alias	AMHR; MRII; MISR2; MISRII; AMHR2
Species	Human
GeneID (Human)	269
GeneID	269
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-144 of human AMHR2 (NP_065434.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q16671
Protein Weight	63 kDa
Shipping	Ice bag