


Anti-Amot Polyclonal Antibody For Wb, Elisa - Q4Vcs5

Item Number: CM0020032



Introduction

Rabbit polyclonal antibody targeting human angiotensin (AMOT), validated for WB and ELISA with cross-reactivity to mouse and rat. Ideal for tight junction and cell migration studies.

[Learn More](#)

Parameter	Value
Product Name	AMOT Rabbit pAb
Remarks/Alias	AMOT; angiotensin
Species	Human
GeneID (Human)	154796
GeneID	154796
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1000-1084 of human AMOT (NP_001106962.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q4VCS5
Protein Weight	118kDa
Shipping	Ice bag