


Anti-Amot Rabbit Polyclonal Antibody For Wb, Elisa - Q4Vcs5

Item Number: CM0025415



Introduction

Rabbit polyclonal antibody targeting human AMOT (angiomotin) for Western blot and ELISA. Suitable for tight junction, endothelial cell migration, and Hippo signaling research. Cross-reacts with rat. SWISS: Q4VCS5.

[Learn More](#)

Parameter	Value
Product Name	AMOT Rabbit pAb
Remarks/Alias	AMOT
Species	Human
GeneID (Human)	154796
GeneID	154796
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human AMOT (NP_573572.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	Q4VCS5
Protein Weight	118kDa
Shipping	Ice bag