

Anti-Amotl1 Rabbit Polyclonal Antibody For Wb And Elisa - Q8Iy63

Item Number: CM0024599

Introduction

Rabbit polyclonal antibody against human AMOTL1 (Angiomotin-like protein 1) validated for Western blot and ELISA. Cross-reacts with mouse and rat. Ideal for Wnt signaling and tight junction research.

[Learn More](#)

Parameter	Value
Product Name	AMOTL1 Rabbit pAb
Remarks/Alias	JEAP; AMOTL1
Species	Human
Gene ID (Human)	154810
Gene ID	154810
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 810-900 of human AMOTL1 (NP_570899.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8IY63
Protein Weight	107kDa
Shipping	Ice bag