

Anti-Amotl2 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Y2J4

Item Number: CM0021677

Introduction

Rabbit polyclonal antibody against AMOTL2, suitable for WB, IHC-P, IF/ICC, ELISA. Detects human and mouse angiomin-like protein 2, a regulator of actin architecture, cell junctions, and Hippo/Wnt signaling.

[Learn More](#)

Parameter	Value
Product Name	AMOTL2 Rabbit pAb
Remarks/Alias	LCCP; AMOTL2
Species	Human
Gene ID (Human)	51421
Gene ID	51421
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 638-837 of human AMOTL2 (NP_001265612.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9Y2J4
Protein Weight	86kDa
Shipping	Ice bag