

Anti-Ampd1 Polyclonal Antibody For Wb, If/Icc, Elisa - P23109

Item Number: CM0024972

Introduction

Rabbit polyclonal antibody against human AMPD1 (AMP deaminase 1), validated in WB, IF/ICC, and ELISA. Suitable for human, mouse, and rat samples; assesses energy metabolism enzyme.

[Learn More](#)

Parameter	Value
Product Name	AMPD1 Rabbit pAb
Remarks/Alias	MAD; MADA; MMDD; AMPD1
Species	Human
GeneID (human)	270
GeneID	270
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 50-260 of human AMPD1 (NP_001166097.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23109
Protein Weight	86kDa
Shipping	Ice bag