


Anti-Amph Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P49418

Item Number: CM0029157



Introduction

Anti-AMPH rabbit polyclonal antibody validated for Western blot, IF/ICC, and ELISA. Detects amphiphysin (~76 kDa) in human, mouse, and rat samples. Ideal for exocytosis, synaptic vesicle, and cytoskeleton research.

[Learn More](#)

Parameter	Value
Product Name	AMPH Rabbit pAb
Remarks/Alias	AMPH1; AMPH
Species	Human
GeneID (Human)	273
GeneID	273
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49418
Protein Weight	76kDa
Shipping	Ice bag