


Anti-Amphiregulin (Areg) Rabbit Polyclonal Antibody For If/Icc, Elisa - P15514

Item Number: CM0014453



Introduction

Rabbit polyclonal antibody targeting human Amphiregulin (AREG), validated for IF/ICC and ELISA. Cross-reacts with mouse. Target weight ~28 kDa. Ideal for studying EGFR ligand signaling and cell growth.

[Learn More](#)

Parameter	Value
Product Name	AREG Rabbit pAb
Remarks/Alias	AR; SDGF; AREGB; CRDGF; AREG
Species	Human
GeneID (Human)	374
GeneID	374
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-252 of human AREG (NP_001648.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Mouse
SWISS	P15514
Protein Weight	28kDa
Shipping	Ice bag