


Anti-Amphiregulin Monoclonal Antibody For Wb, Elisa - P15514

Item Number: CM0011546



Introduction

Recombinant human Amphiregulin monoclonal antibody validated for WB and ELISA.

Recognizes human and rat Amphiregulin (AREG), a ligand of the EGFR receptor. Ideal for growth factor research.

[Learn More](#)

Parameter	Value
Product Name	Amphiregulin Mouse mAb
Remarks/Alias	AR; SDGF; AREGB; CRDGF
Species	Human
GeneID (Human)	374
GeneID	374
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 25-252 of human Amphiregulin (NP_001648.1).
Source	Mouse
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P15514
Protein Weight	28 kDa
Shipping	Ice bag