

Anti-Ampk Alpha 1/2 Monoclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - Q13131 / P54646

Item Number: CM0012109



Introduction

Rabbit monoclonal antibody against human AMPK catalytic subunits alpha-1 (PRKAA1) and alpha-2 (PRKAA2). Validated for WB, IF/ICC, IHC-P, and ELISA, cross-reacts with mouse and rat. Suitable for metabolic, autophagy, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	AMPK alpha 1/2 Rabbit mAb
Remarks/Alias	AMPK; AMPKa1; AMPK alpha 1; AMPK alpha 2; AMPKa2;
Species	Human
Gene ID (Human)	5562/5563
Gene ID	5562/5563
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 13-300 of human AMPK alpha 1 (NP_006242.5).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13131/P54646
Protein Weight	64kDa
Shipping	Ice bag