

Anti-Ampk α 1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q13131

Item Number: CM0016130



Introduction

Rabbit polyclonal antibody against AMPK α 1 (PRKAA1), validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Detects endogenous AMPK α 1 protein at 64 kDa.

[Learn More](#)

Parameter	Value
Product Name	AMPK α 1 Rabbit pAb
Remarks/Alias	AMPK; AMPK α 1; AMPK alpha 1; AMPK α 1
Species	Human
GeneID	5562
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 381-520 of human AMPK α 1 (NP_006242.5).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13131
Protein Weight	64kDa
Shipping	Ice bag