

Anti-Ampkα1 Rabbit Polyclonal Antibody For Wb, Elisa - Q13131

Item Number: CM0019051



Introduction

High-quality rabbit polyclonal antibody targeting AMPKα1 (PRKAA1), validated for Western blot and ELISA with mouse cross-reactivity. AMPKα1 is a central energy sensor regulating metabolism, autophagy, and transcription—ideal for research and diagnostic applications.

[Learn More](#)

Parameter	Value
Product Name	AMPKα1 Rabbit pAb
Remarks/Alias	AMPK; AMPKα1; AMPK alpha 1; AMPKα1
Gene ID (Human)	5562
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-200 of human AMPKα1 (NP_006242.5).
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q13131
Protein Weight	64kDa
Shipping	Ice bag