

Anti-Ampka1/Ampka2 Monoclonal Antibody For Wb, Ihc-P, Elisa - Q13131 / P54646

Item Number: CM0012014



Introduction

Mouse monoclonal antibody against human AMPK α 1/ α 2 (Q13131/P54646). Validated for WB, IHC-P, ELISA. Cross-reacts with Mouse, Rat. Ideal for energy metabolism, autophagy, and cancer signaling research.

[Learn More](#)

Parameter	Value
Product Name	AMPK α 1/AMPK α 2 Mouse mAb
Remarks/Alias	AMPK α 1/AMPK α 2
Species	Human
Gene ID (Human)	5562/5563
Gene ID	5562/5563
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 251-550 of human AMPK α 1 (Q13131)
Source	Mouse
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13131/P54646
Protein Weight	64kDa/65kDa
Shipping	Ice bag